Floorline LTC











Our Floorline LTC is the optimal bed solution for home care use or long-term care facilities. Its one piece design makes it easy to fold, store and transport. The Floorline LTC comes in two convenient options including one with a wood-grain finish that compliments any room.

Clinical and Patient Benefits

The Floorline LTC lowers to a distance of 9.5cm/3.7" from the floor. This distance, although seemingly small, has enormous benefits to both the patient and the facility. By just adding 5cm to the distance between the bed and the floor, the impact force of a fall is increased by 50%. This can be the difference between a simple bruise or a fractured hip for the patient that may result in an extented hospital stay or even death. In addition to the physical effects, fall injuries also have psychological effects on the patient as they experience a loss of confidence and fear after a fall.

Optimal Features

The Floorline LTC features standard positions

normally found only on hospital beds, such as Trendelenburg, reverse Trendelenburg and cardiac positioning. With a full assortment of accessories, an easy-to-use handset, and durable LINAK electronics, the Floorline LTC is the best choice for patient care.

Quality

Our products are designed and manufactured to meet all the requirements of relevant standards. At every stage of the development process, every single product is given extensive attention in design, development and manufacturing. This means that a user of this product feels confident with its reliability and the caregiver feels confident with its safety.











Accessories



IV Pole Adjustable height. Can be mounted on either end of bed. Maximum SWL 7kg/ 15.4lb. Part no. IVP



Care Assist Rail Designed to assist patients safely in and out of a bed. Ergonomic design. Part no. 90702



Extension Kit Extends bed by 17.5cm/7" for taller patients. Part no. 90709



Floorline LTC

Handset Five function - included with bed. Part no. 75746 Handset holder Part no. LC2250, Lock Out Key Part no. 75790

Technical Specifications

, interest	()
1	I
	1
1	





at 1

Wall Bumper Bar Protects the head of the bed. Mounted on lower support bracket between wheels. Part no. 90705

Self Help Pole

SHP600

Assists patient moving

within confines of

bed. Maximum SWL

75kg/165.3lb. Part no.



Oxygen Holder Carries C size tank. Attaches to legs. Part no.





Compact Folds compactly for

Standard Positions backrest, kneebreak, Trendelenburg, reverse Trendelenburg, cardiac.



Premium Casters One central locking and three steering/directional casters for premium

maneuvrability.

Includes auto-contour,





Part no. White Finish 90600 (EU), 90602 (UK), 90601 (EU), 90603 (UK), Part no. Woodgrain Finish 90606 (NA), 90608 (AU) 90607(NA), 90609(AU) Platform Position 9.5cm - 70cm / 3.7" - 27.5" 215cm - 232.5cm / Overall Bed Length (Standard to (Min-Max) Extended with Extension Kit) 84.5" - 91.5" Platform Dimensions (LxW) 215cm x 90.5cm / Platform Dimensions with 232.5cm x 90.5cm / 84.5" x 35.5" Extension Kit $(L \times W)$ 91.5" x 35.5" Sleeping Dimensions (LxW) 200cm x 90cm / Sleeping Dimensions with 217.5cm x 90cm / 78.7" x 35.5" Extension Kit (LxW) 85.5" x 35.5" Output Voltage Max. 24V DC IP Rating IPX4 Input Voltage 120V, 60Hz (90404 to 90607), 106kg/234lb Product Weight else 230V, 50Hz Platform Panel Angles Backrest 70°, Thigh 45°, Knee Safe Working Load (SWL) 200kg/440lb 110°, Calf 25°, Trendelenburg/ Reverse Trendelenburg 18° Duty Cycle 2 minutes operation, 18 minutes Maximum User Weight 165kg/363lb rest